




LDB3 Antibody

Product Code	CSB-PA012831ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O75112
Immunogen	Recombinant Human LIM domain-binding protein 3 protein (1-283AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM domain-binding protein 3 (Protein cypher) (Z-band alternatively spliced PDZ-motif protein), LDB3, KIAA0613 ZASP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LDB3

Image



Western blot

All lanes: LDB3 antibody at 3.21µg/ml

Lane 1 : Hela whole cell lysate

Lane 2 : Raji whole cell lysate

Lane 3 : MCF-7 whole cell lysate

Lane 4 : A431 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 78, 67, 51, 43, 36, 31, 79

kDa

Observed band size: 67, 51 kDa



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA012831ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.