





ANKS1B Antibody

Product Code	CSB-PA767207ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z6G8
Immunogen	Recombinant Human Ankyrin repeat and sterile alpha motif domain-containing protein 1B protein (999-1248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ankyrin repeat and sterile alpha motif domain-containing protein 1B (Amyloid-beta protein intracellular domain-associated protein 1) (AIDA-1) (E2A-PBX1-associated protein) (EB-1), ANKS1B
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ANKS1B

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA767207ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA767207ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.