





CACNG7 Antibody

Product Code	CSB-PA004421ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62955
Immunogen	Recombinant Human Voltage-dependent calcium channel gamma-7 subunit protein (156-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent calcium channel gamma-7 subunit (Neuronal voltage-gated calcium channel gamma-7 subunit) (Transmembrane AMPAR regulatory protein gamma-7) (TARP gamma-7), CACNG7
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CACNG7

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA004421ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA004421ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.