




CHRNG Antibody

Product Code	CSB-PA005401ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07510
Immunogen	Recombinant Human Acetylcholine receptor subunit gamma protein (328-517AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acetylcholine receptor subunit gamma, CHRNG, ACHRG
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CHRNG

Image



Western blot

All lanes: CHRNG antibody at 5.5µg/ml

Lane 1: 293T whole cell lysate

Lane 2: HepG2 whole cell lysate

Lane 3: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 58, 53 kDa

Observed band size: 58 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.