





FGF12 Antibody

Product Code	CSB-PA008618ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61328
Immunogen	Recombinant Human Fibroblast growth factor 12 protein (1-181AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 12 (FGF-12) (Fibroblast growth factor homologous factor 1) (FHF-1) (Myocyte-activating factor), FGF12, FGF12B FHF1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGF12

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA008618ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA008618ESR1HU at dilution of 1:100