




VSX1 Antibody

Product Code	CSB-PA025938ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZR4
Immunogen	Recombinant Human Visual system homeobox 1 protein (1-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Visual system homeobox 1 (Homeodomain protein RINX) (Retinal inner nuclear layer homeobox protein) (Transcription factor VSX1), VSX1, RINX
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VSX1

Image



Western blot

All lanes: VSX1 antibody IgG at 2.32µg/ml + 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 39, 25, 40, 15, 30, 23, 25, 33 kDa

Observed band size: 39 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.