





GBF1 Antibody

Product Code	CSB-PA846077ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92538
Immunogen	Recombinant Human Golgi-specific brefeldin A-resistance guanine nucleotide exchange factor 1 protein (1-85AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Golgi-specific brefeldin A-resistance guanine nucleotide exchange factor 1 (BFA-resistant GEF 1), GBF1, KIAA0248
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GBF1


Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA846077ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA846077ESR1HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA846077ESR1HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.