




ARHGEF3 Antibody

Product Code	CSB-PA878863ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR81
Immunogen	Recombinant Human Rho guanine nucleotide exchange factor 3 protein (287-526AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho guanine nucleotide exchange factor 3 (Exchange factor found in platelets and leukemic and neuronal tissues) (XPLN), ARHGEF3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ARHGEF3

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA878863ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.