





PPIL2 Antibody

Product Code	CSB-PA621661ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13356
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase-like 2 protein (1-260AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase-like 2 (PPIase) (EC 2.3.2.27) (EC 5.2.1.8) (CYC4) (Cyclophilin-60) (Cyclophilin-like protein Cyp-60) (Cyp60) (hCyp-60) (RING-type E3 ubiquitin transferase isomerase-like 2) (Rotamase PPIL2), PPIL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPIL2

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA621661ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA621661ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.