




UBE2F Antibody

Product Code	CSB-PA850248ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969M7
Immunogen	Recombinant Human NEDD8-conjugating enzyme UBE2F protein (1-185AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NEDD8-conjugating enzyme UBE2F (EC 2.3.2.-) (NEDD8 carrier protein UBE2F) (NEDD8 protein ligase UBE2F) (NEDD8-conjugating enzyme 2) (RING-type E3 NEDD8 transferase UBE2F) (Ubiquitin-conjugating enzyme E2 F), UBE2F, NCE2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UBE2F

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA850248ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.