



MESP2 Antibody

Product Code	CSB-PA613150ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0VG99
Immunogen	Recombinant Human Mesoderm posterior protein 2 protein (218-397AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mesoderm posterior protein 2 (Class C basic helix-loop-helix protein 6) (bHLHc6), MESP2, BHLHC6 SCDO2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MESP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.