





GATC Antibody

Product Code	CSB-PA009283ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43716
Immunogen	Recombinant Human Glutamyl-tRNA(Gln) amidotransferase subunit C, mitochondrial protein (1-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamyl-tRNA(Gln) amidotransferase subunit C, mitochondrial (Glu-AdT subunit C) (EC 6.3.5.-) (Protein 15E1.2), GATC, 15E1.2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GATC

Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA009283ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA009283ESR1HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.