




NFKBIE Antibody

Product Code	CSB-PA015764ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00221
Immunogen	Recombinant Human NF-kappa-B inhibitor epsilon protein (140-280AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NF-kappa-B inhibitor epsilon (NF-kappa-BIE) (I-kappa-B-epsilon) (IκB-E) (IκB-epsilon) (IkappaBepsilon), NFKBIE, IKBE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NFKBIE

Image



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA015764ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.