




ATXN10 Antibody

Product Code	CSB-PA883357ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBB4
Immunogen	Recombinant Human Ataxin-10 protein (346-475AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ataxin-10 (Brain protein E46 homolog) (Spinocerebellar ataxia type 10 protein), ATXN10, SCA10
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ATXN10

Image



Western blot

All lanes: ATXN10 antibody at 4.79μg/ml

Lane 1: A431 whole cell lysate

Lane 2: LO2 whole cell lysate

Lane 3: HepG2 whole cell lysate

Lane 4: 293T whole cell lysate


Lane 5: Hela whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 54, 47 kDa


Observed band size: 54 kDa



Immunohistochemistry of paraffin-embedded

human testis tissue using CSB-

PA883357ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA883357ESR1HU at dilution of 1:100