




GSTK1 Antibody

Product Code	CSB-PA009978ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2Q3
Immunogen	Recombinant Human Glutathione S-transferase kappa 1 protein (1-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase kappa 1 (EC 2.5.1.18) (GST 13-13) (GST class-kappa) (GSTK1-1) (hGSTK1) (Glutathione S-transferase subunit 13), GSTK1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GSTK1

Image



Western blot

All lanes: GSTK1 antibody at 0.51µg/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: K562 whole cell lysate

Lane 3: Hela whole cell lysate

Lane 4: HepG2 whole cell lysate


Lane 5: Raji whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 26, 32, 25, 21 kDa


Observed band size: 26 kDa



Immunohistochemistry of paraffin-embedded

human liver cancer using CSB-

PA009978ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA009978ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.