




NSMF Antibody

Product Code	CSB-PA859576ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6X4W1
Immunogen	Recombinant Human NMDA receptor synaptonuclear signaling and neuronal migration factor protein (251-530AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	NMDA receptor synaptonuclear signaling and neuronal migration factor (Nasal embryonic luteinizing hormone-releasing hormone factor) (Nasal embryonic LHRH factor), NSMF, NELF
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NSMF

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA859576ESR1HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.