




PPIF Antibody

Product Code	CSB-PA018476ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30405
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase F, mitochondrial protein (30-207AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase F, mitochondrial (PPIase F) (EC 5.2.1.8) (Cyclophilin D) (CyP-D) (CypD) (Cyclophilin F) (Mitochondrial cyclophilin) (CyP-M) (Rotamase F), PPIF, CYP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PPIF

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA018476ESR2HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.