




SH3GL3 Antibody

Product Code	CSB-PA860359ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99963
Immunogen	Recombinant Human Endophilin-A3 protein (1-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endophilin-A3 (EEN-B2) (Endophilin-3) (SH3 domain protein 2C) (SH3 domain-containing GRB2-like protein 3), SH3GL3, CNSA3 SH3D2C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SH3GL3

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA860359ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.