




GORASP2 Antibody

Product Code	CSB-PA880999DSR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8Y8
Immunogen	Recombinant Human Golgi reassembly-stacking protein 2 protein (353-452AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Golgi reassembly-stacking protein 2 (GRS2) (Golgi phosphoprotein 6) (GOLPH6) (Golgi reassembly-stacking protein of 55 kDa) (GRASP55) (p59), GORASP2, GOLPH6
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GORASP2

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA880999DSR2HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.