




STC2 Antibody

Product Code	CSB-PA022822ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76061
Immunogen	Recombinant Human Stanniocalcin-2 protein (25-302AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Stanniocalcin-2 (STC-2) (Stanniocalcin-related protein) (STC-related protein) (STCRP), STC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	STC2

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA022822ESR1HU at dilution of 1:100