




TMEFF2 Antibody

Product Code	CSB-PA883439ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIK5
Immunogen	Recombinant Human Tomoregulin-2 protein (41-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tomoregulin-2 (TR-2) (Hyperplastic polyposis protein 1) (Transmembrane protein with EGF-like and two follistatin-like domains), TMEFF2, HPP1 TENB2 TPEF
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	TMEFF2

Image



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA883439ESR2HU at dilution of 1:100