




P2RX3 Antibody

Product Code	CSB-PA017321ESR1HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P56373
Immunogen	Recombinant Human P2X purinoceptor 3 protein (44-317AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	P2X purinoceptor 3 (P2X3) (ATP receptor) (Purinergic receptor), P2RX3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	P2RX3

Image



Western blot


All lanes: P2RX3 antibody at 3.91 µg/ml +
Mouse liver tissue

Secondary


Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 44 kDa

Observed band size: 44 kDa



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA017321ESR1HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human pancreatic tissue using CSB-PA017321ESR1HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.