



mcr1 Antibody, HRP conjugated

Product Code	CSB-PA745804LB01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A0R6L508
Immunogen	Recombinant Escherichia coli Probable phosphatidylethanolamine transferase Mcr-1 protein (1-541AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable phosphatidylethanolamine transferase Mcr-1 (EC 2.7.-.-) (Polymyxin resistance protein MCR-1), mcr1, mcr-1
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	mcr1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.