





ULK4 Antibody

Product Code	CSB-PA853394ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C45
Immunogen	Recombinant Human Serine/threonine-protein kinase ULK4 protein (1-290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase ULK4 (EC 2.7.11.1) (Unc-51-like kinase 4), ULK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ULK4

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA853394ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human skin tissue using CSB-PA853394ESR2HU at dilution of 1:100