




Recombinant Human Legumain (LGMN)

Product Code	CSB-EP012903HU
Relevance	Has a strict specificity for hydrolysis of asparaginyl bonds. Can also cleave aspartyl bonds slowly, especially under acidic conditions. Required for normal lysosomal protein degradation in renal proximal tubules. Required for normal degradation of internalized EGFR. Plays a role in the regulation of cell proliferation via its role in EGFR degradation (By similarity). May be involved in the processing of proteins for MHC class II antigen presentation in the lysosomal/endosomal system.
Abbreviation	Recombinant Human LGMN protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q99538
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	VPIDDPEDGGKHWVIVAGSNGWYNYRHHQADACHAYQIIHRNGIPDEQIVVM MYDDIAYSEDNPTPGIVINRPNNGTDVYQQGVPKDYTGEDVTPQNFLAVLRGDAE AVKGIGSGKVLKSGPQDHFVFIYFTDHGSTGILVFPNEDLHVKDLNETIHMYMKH KMYRKMVFYIEACESGSMNHLDPNINVYATTAANPRESSYACYDEKRSTYL GDWYSVNW MEDSDVEDLTKETLHKQYHLVKSHNTNTSHVMQYGNKTISTMKV MQFQGMKRKASSPVPLPPVTHLDLTPSPDVPLTIMKRKLMNTN
Research Area	Immunology
Source	E.coli
Target Names	LGMN
Protein Names	Recommended name: Legumain EC= 3.4.22.34Alternative name(s): Asparaginyl endopeptidase Protease, cysteine 1
Expression Region	18-323aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	38.8kDa
Protein Length	Full Length of Mature Protein of BC003061
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.