




Recombinant Mouse Zinc finger protein 346 (Znf346)

Product Code	CSB-EP026691MO
Relevance	Binds with low affinity to dsDNA and ssRNA, and with high affinity to dsRNA, with no detectable sequence specificity.
Abbreviation	Recombinant Mouse Znf346 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9R0B7
Alias	Just another zinc finger protein
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MECPAPDATDAADPGEAGPYKGSEEPEGREPDGVRFDREARRRLWEAVSGA QPAGREEVEHMIQKNQCLFTSTQCKVCCAMLISESQKLAHYQSKKHANKVKR YLAIHGMETIKGDVKRLSDQKSSRSKDKNHCCPICNMTFSSPAVAQSHYLGK THAKSLKLKQQSTKGAALQQNREMLDPDKFCSLCHSTFNDPAMAQQHYMGK RHRKQETKLKLMAHYGRLADPAVSDLPAGKGYPCCKTCKIVLNSIEQYQAHVSG FKHKNQSPKTLVTLGSQTPVQTQPTPKDSSTVQD
Research Area	Others
Source	E.coli
Target Names	Znf346
Expression Region	1-294aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	48.7kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.