




Recombinant Human PHD finger protein 11 (PHF11)

Product Code	CSB-EP017900HU
Relevance	Positive regulator of Th1-type cytokine gene expression.
Abbreviation	Recombinant Human PHF11 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UIL8
Alias	BRCA1 C-terminus-associated protein Renal carcinoma antigen NY-REN-34
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MEKRTCALCPKDVEYNVLYFAQSENIAAHENCLLYSSGLVECEDQDPLNPDRS FDVESVKKEIQRGRKLKCKFCHKRGATVGCCLKNCNKNYHFFCAKKDDAVPQ SDGVRGIYKLLCQQHAQFPPIAQSASFSGVKRKRGRKKPLSGNHVQPPETMKC NTFIRQVKKEEHGRHTDATVKVPFLKKCKEAGLLNYLLEEILDKVHSIPEKLMDDET TSESDYEEIGSALFDCRLFEDTFVNFQAAIEKKIHASQQRWQQLKEEIELLQDL KQTLCSFQENRDLMSSTSISSLSY
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	PHF11
Protein Names	Recommended name: PHD finger protein 11Alternative name(s): BRCA1 C-terminus-associated protein Renal carcinoma antigen NY-REN-34
Expression Region	1-292aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	60.5kDa
Protein Length	Full Length of Isoform 2
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.