




Recombinant Mouse Tenomodulin (Tnmd), partial

Product Code	CSB-EP024007MO
Abbreviation	Recombinant Mouse Tnmd protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9EP64
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KHFWPEVSKKTYDMEHTFYSNGEKKKIYMEIDPITRTEIFRSGNGTDETLVHD FKNGYTGIVFVGLQKCFIKTQIKVIPEFSEPEEEIDENEEITTTFFEQSVIWP KPIENRDFLKNKILEICDNVTMYWINPTLIAVSELQDFEEDGEDLHFPTSEKKGI DQNEQWVVPQVKVEKTRHTRQASEEDLPINDYTENGIEFDPMLDERGYCCY CRRGNRYCRRVCEPLLGYPPYPYCYQGGRVICRVIMPCNWWWARMMLGRV
Research Area	Others
Source	E.coli
Target Names	Tnmd
Protein Names	Recommended name: Tenomodulin Short name= TeM Short name= mTeM Alternative name(s): Chondromodulin-1-like protein Short name= ChM1L Short name= mChM1L Chondromodulin-I-like protein Myodulin Tendin
Expression Region	51-317aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	35.6 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.