




Recombinant Human Zinc finger protein 155 (ZNF155)

Product Code	CSB-EP026546HU
Relevance	May be involved in transcriptional regulation.
Abbreviation	Recombinant Human ZNF155 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12901
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MTTFKEAVTFKDVAVVFTEEEELGLLDPAQRKLYRDVMLENFRNLLSVGHQPFH QDTCHFLREEKFWMMGTATQREGNSGGKIQTELESVPEAGAHEEWSCQQIW EQIAKDLTRSQDSIINNSQFFENGDVPSQVEAGLPTIHTGQKPSQGGKCKQSSIS DVPIFDLPQQLYSEEKSYTCDECGKSICYISALHVVHQRVHVGEKLFMCDVCGK EFSQSSHLQTHQRVHTGEKPFKCEQCGKGFSRRSALNVHRKLHTGEKPYICE ACGKAIFIHDSQLKEHKRIHTGEKPFKCDICGKTFYFRSRLKSHSMVHTGEKPF RCDTCDKSFHQRSALNRHCMVHTGEKPYRCEQCGKGFIGRLDFYKHQVVHT GEKPYNCKECKGKSFWRSSCLLNHQRVHSGEKSFKCEECKGKGYTNSQLSSH QRSHSGEKPYKCEECKGKYVTKFNLDLHQRVHTGERPYNCKECKGNFSRAS SILNHLKRLHCQKKPFKCEDCGKRLVHRTYRKDQPRDYSGENPSKCEDCGRRY KRRLNLDILLSLFLNDT
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	ZNF155
Expression Region	1-538aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	78.1kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.