




Recombinant Dog Annexin A13 (ANXA13)

Product Code	CSB-MP001839DO
Relevance	Involved in vesicular traffic to the apical plasma membrane (isoform B).
Abbreviation	Recombinant Dog ANXA13 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q29471
Product Type	Recombinant Proteins
Immunogen Species	Canis lupus familiaris (Dog)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GNRHAKAKSHHGFDVDHDAKLNKACKGMGTDEAAIIEILSSRTSDERQQIKQ KYKATYGKDLEEVFKSDLSGNFEKTALALLDRPSEYDARQLQKAMKGLGTDEA VLIEILCTR TNKEIMAIKEAYQRLFDRSLESDVKADTSGNLKAILVSLQANRDEG DDVDKDLAQDAKDLYDAGDGRWGTDELAFNEVLAKRSHKQLRATFQAYQIL IDKDIEEAIEAETSGDLQKAYLTLVRCARDQEGYFADRLYKSMKGTGTDEETLI HIIVTRAEVDLQGIKAKFQEKYQKSLSDMVRSDTSGDFQKLLVALLH
Research Area	Signal Transduction
Source	Mammalian cell
Target Names	ANXA13
Protein Names	Recommended name: Annexin A13Alternative name(s): Annexin XIII Annexin-13 Intestine-specific annexin Short name= ISA
Expression Region	2-316aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Myc-tagged
Mol. Weight	39.3kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.