




Recombinant *Macaca fascicularis* Legumain (LGMN)

Product Code	CSB-EP012903MOV
Relevance	Has a strict specificity for hydrolysis of asparaginyl bonds. Can also cleave aspartyl bonds slowly, especially under acidic conditions. Required for normal lysosomal protein degradation in renal proximal tubules. Required for normal degradation of internalized EGFR. Plays a role in the regulation of cell proliferation via its role in EGFR degradation. May be involved in the processing of proteins for MHC class II antigen presentation in the lysosomal/endosomal system
Abbreviation	Recombinant Cynomolgus monkey LGMN protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4R4T8
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (Cynomolgus monkey)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	VPIDDPEDGGKHWVIVAGSNGWYNYRHQADACHAYQIIHRNGIPDEQIVVM MYDDIAYSEDNPTPGIVINRPNGTDVYQGVPKDYTGEDVTPQNFLAVLRGDAE AVKGIGSGKVLKSGPQDHFVYFTDHGSGTILVFPNEDLHVKDLNETIYYMYKH KMYRKMVFYIEACESGSMNHLPDNINVYATTAANPRESSYACYDEKRSTYL GDWYSVNW MEDSDVEDLTKETLHKQYHLVKSHTNTSHVMQYGNKTIISTMKV MQFQGMKHKASSPLSLPPVTHLDLTPSPDVPLTIMKRKLMNTN
Research Area	Immunology
Source	<i>E. coli</i>
Target Names	LGMN
Protein Names	Asparaginyl endopeptidase Protease, cysteine 1 Gene names: LGMN PRSC1
Expression Region	18-323aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	39.7 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.