



Recombinant Human Metalloproteinase inhibitor 3 (TIMP3)

Product Code	CSB-EP023562HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35625
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	CTCSPSH PQDAFCNSDI VIRAKVVGKK LVKEGPFGLT VYTIKQMKMY RGFTKMPHVQ YIHTEASESL CGLKLEVNKY QYLLTGRVYD GKMYTGLCNF VERWDQLTLS QRKGLNYRYH LGCNCKIKSC YYLPCFVTSK NECLWTDMLS NFGYPGYQSK HYACIRQKGG YCSWYRGWAP PDKSIINATD P
Source	E.coli
Target Names	TIMP3
Protein Names	Recommended name: Metalloproteinase inhibitor 3 Alternative name(s): Protein MIG-5 Tissue inhibitor of metalloproteinases 3 Short name= TIMP-3
Expression Region	24-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene belongs to the TIMP gene family. The proteins encoded by this gene family are inhibitors of the matrix metalloproteinases, a group of peptidases involved in degradation of the extracellular matrix (ECM). Expression of this gene is induced in response to mitogenic stimulation and this netrin domain-containing protein is localized to the ECM. Mutations in this gene have been associated with the autosomal dominant disorder Sorsby s fundus dystrophy.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.