



# Recombinant *Drosophila erecta* Tryptophan 2,3-dioxygenase (v)

<b>Product Code</b>	CSB-EP023351DLK-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	B3NVC6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Drosophila erecta</i> (Fruit fly)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>                     MSCPYAGNGN DHDDSAVPLT TEVGKIYGEY LMLDKLLDAQ CMLSEEDKRP                      VHDEHLFIIT HQAYELWFKQ IIFEFDISIRD MLDAEVIDET KTLEIVKRLN                      RVVLILKLLV DQVPILETMT PLDFMDFRKY LAPASGFQSL QFRLIENKLG                      VLTEQVRVRYN QKYSDFVFSDE EARNSIRNSE KDPSLLELVQ RWLERTPGLE                      ESGFNFWAKF QESVDRFLEA QVQSAMMEPV EKAKNYRLMD IEKRREYVRS                      IFDPAVHDAL VRRGDRRFSH RALQGAIMIT FYRDEPRFSQ PHQLLTLLMD                      IDSLITKWRY NHVIMVQRMI GSQQLGTGGS SGYQYLRSTL SDRYKVFLDL                      FNLSTFLIPR EAIPPLDETI RRKLVHKS V                 </p>
<b>Source</b>	<i>E.coli</i>
<b>Target Names</b>	v
<b>Protein Names</b>	Recommended name: Tryptophan 2,3-dioxygenase Short name= TDO EC= 1.13.11.11 Alternative name(s): Protein vermilion Tryptamin 2,3-dioxygenase Tryptophan oxygenase Short name= TO Short name= TRPO Tryptophan py
<b>Expression Region</b>	1-379
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.