



Recombinant *Drosophila erecta* Tryptophan 2,3-dioxygenase (v)

Product Code	CSB-EP023351DLK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B3NVC6
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila erecta</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	<p> MSCPYAGNGN DHDDSAVPLT TEVGKIYGEY LMLDKLLDAQ CMLSEEDKRP VHDEHLFIIT HQAYELWFKQ IIFEFDISIRD MLDAEVIDET KTLEIVKRLN RVVLILKLLV DQVPILETMT PLDFMDFRKY LAPASGFQSL QFRLIENKLG VLTEQRVRYN QKYSDFVFSDE EARNSIRNSE KDPSLLELVQ RWLERTPGLE ESGFNFWAKF QESVDRFLEA QVQSAMMEPV EKAKNYRLMD IEKRREYVRS IFDPAVHDAL VRRGDRRFSSH RALQGAIMIT FYRDEPRFSQ PHQLLTLLMD IDSLITKWRY NHVIMVQRMI GSQQLGTGGS SGYQYLRSTL SDRYKVFLDL FNLSTFLIPR EAIPPLDETI RRKLVHKS V </p>
Source	<i>E.coli</i>
Target Names	v
Protein Names	Recommended name: Tryptophan 2,3-dioxygenase Short name= TDO EC= 1.13.11.11 Alternative name(s): Protein vermilion Tryptamin 2,3-dioxygenase Tryptophan oxygenase Short name= TO Short name= TRPO Tryptophan py
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.