



Recombinant Escherichia coli O157:H7 Argininosuccinate lyase (argH)

Product Code	CSB-EP002213EOE-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B5Z062
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7 (strain EC4115 / EHEC)
Purity	>85% (SDS-PAGE)
Sequence	MALWGGRFTQ AADQRFKQFN DSLRFDYRLA EQDIVGSVAW SKALVTVGVL TAEFAQLEE ALNVLLEDVR ARPQQILESD AEDIHSWVEG KLIDKVGQLG KKLHTGRSRN DQVATDLKLW CKDTVSELLT ANRQLQSALV ETAQNNQDAV MPGYTHLQRA QPVTFAHWCL AYVEMLARDE SRLQDALKRL DVSPLGCGAL AGTAYEIDRE QLAGWLG FAS ATRNSLDSVS DRDHVLELLS AAAIGMVHLS RFAEDLIFFN TGEAGFVELS DRVTS GSSLM PQKKNPDALE LIRGKCGRVQ GALTGMMMTL KGLPLAYNKD MQEDKEGLFD ALDTWLDCLH MAALVLDGIQ VKRPRCQEAQ QQGYANATEL ADYLVAKGVP FREAHHIVGE AVVEAIRQ GK PLEELPLTEL QKFSPVIGED VYPILSLQSC LDKRAAKGGV SPQQVAQAIA FARARLG
Source	E.coli
Target Names	argH
Protein Names	Recommended name: Argininosuccinate lyase Short name= ASAL EC= 4.3.2.1 Alternative name(s): Arginosuccinase
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.