



Recombinant Influenza A virus Non-structural protein 1 (NS)

Product Code	CSB-MP604972IER
Abbreviation	NS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0HD54
Product Type	Recombinant Protein
Immunogen Species	Influenza A virus (strain A/Hickox/1940 H1N1)
Purity	>85% (SDS-PAGE)
Sequence	MDPNTVSSFQ VDCFLWHVRK RVADQELGDA PFLDRLRRDQ KSLKGRGSTL GLNIETATRV GKQIVERILK EESDEALKMT MASAPASRYL TDMTIEEMSR DWFMLMPKQK VAGPLCIRMD QAVMDKSIL KANFSVIFDR LETLILLRAF TEEGAIVGEI SPLPSLPGHT NEDVKNAIGV LIGGLEWNDN TVRVSKTLQR FAWRSSNENG GPPLTPKQKR KMARTIRSEV RRNKMVD
Source	Mammalian cell
Target Names	NS
Protein Names	Recommended name: Non-structural protein 1 Short name= NS1 Alternative name(s): NS1A
Expression Region	1-237
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.