



Recombinant Mouse Folate receptor alpha (Folr1)

Product Code	CSB-YP008784MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P35846
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	TRARTE LLNVCMDAKH HKEKPGPEDN LHDQCSPWKT NSCCSTNTSQ EAHKDISYLY RFNWNHCGTM TSECKRHFQ DTCLYECSPN LGPWIQQVDQ SWRKERILDV PLCKEDCQQW WEDCQSSFTC KSNWHKGWNW SSGHNECPVG ASCHPFTFYF PTSAALCEEI WSHSYKLSNY SRGSGRCIQM WFDPAQGNPN EEVARFYAEA MS
Source	Yeast
Target Names	Folr1
Protein Names	Recommended name: Folate receptor alpha Short name= FR-alpha Alternative name(s): Folate receptor 1 Folate-binding protein 1
Expression Region	25-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is a member of the folate receptor family. Members of this gene family bind folic acid and its reduced derivatives, and transport 5-methyltetrahydrofolate into cells. This gene product is a secreted protein that either anchors to membranes via a glycosyl-phosphatidylinositol linkage or exists in a soluble form. Mutations in this gene have been associated with neurodegeneration due to cerebral folate transport deficiency. Due to the presence of two promoters, multiple transcription start sites, and alternative splicing, multiple transcript variants encoding the same protein have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.