




Recombinant Human Killer cell immunoglobulin-like receptor 2DS3 (KIR2DS3), partial

Product Code	CSB-BP613530HU
Relevance	Receptor on natural killer (NK) cells for HLA-C alleles. Does not inhibit the activity of NK cells.
Abbreviation	Recombinant Human KIR2DS3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14952
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	HEGFRRKPSLLAHPGRLVKSEETVILQCWSDVMFEHFLLHREGTFNDTLRLIG EHIDGVSKANFSIGRMRQDLAGTYRCYGSVPHSPYQFSAPSDPLDIVITGLYEK PSLSAQPGPTVLAGESTLSCSSWSSYDMYHLSTEGEAHERFSAAGPKVNGT FQADFPLGPATQGGTYRCFGSFHDSPEYWSKSSDPLLVSVTGNPSNSWPSP TEPSSKTGNPRHLH
Research Area	Others
Source	Baculovirus
Target Names	KIR2DS3
Protein Names	MHC class I NK cell receptor Natural killer-associated transcript 7 Short name: NKAT-7 NKAT7
Expression Region	22-245aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal MBP-tagged and C-terminal 6xHis-tagged
Mol. Weight	68.7 kDa
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.