




# Recombinant *Saccharomyces cerevisiae* Ubiquitin-like protein SMT3 (SMT3)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP615489SVG   |
| <b>Relevance</b>         | Not known; suppressor of MIF2 mutations.  |
| <b>Abbreviation</b>      | Recombinant <i>Saccharomyces cerevisiae</i> SMT3 protein  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| <b>Uniprot No.</b>       | Q12306  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | <i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)  |
| <b>Purity</b>            | ≥ 90% as determined by SDS-PAGE.  |
| <b>Sequence</b>          | SDSEVNQEAKPEVKPEVKPETHINLKVSDGSSEIFFKIKKTTPLRRLMEAFKR<br>QGKEMDSLRLFLYDGIQADQTPEDLDMEDNDIIEAHREQIGG   |
| <b>Research Area</b>     | Others  |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | SMT3  |
| <b>Protein Names</b>     | Recommended name: Ubiquitin-like protein SMT3   |
| <b>Expression Region</b> | 2-98aa  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | N-terminal 6xHis-SUMO-tagged  |
| <b>Mol. Weight</b>       | 27.1kDa   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.