




Recombinant Human Four and a half LIM domains protein 3 (FHL3)

Product Code	CSB-EP618796HU
Abbreviation	Recombinant Human FHL3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13643
Alias	Skeletal muscle LIM-protein 2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SESFDCAKCNESLYGRKYIQTDSGPYCVPCYDNTFANTCAECQQLIGHDSREL FYEDRHFHEGCFRCCRCQRSLADEPFTCQDSELLCND CYCSAFSSQCSACG ETVMPGSRKLEYGGQTWHEHCFLCSGCEQPLGSRSFVPDKGAHYCVPCYEN KFAPRCARCSKTLTQGGVTYRDQPWHRECLVCTGCQTPLAGQQFTSRDEDP YCVACFGELFAPKCSSCKRPIVGLGGGKYVSFEDRHHWHNCFSCARCSTSLV GQGFVPDGDQVLCQGCSQAGP
Research Area	Tags & Cell Markers
Source	E.coli
Target Names	FHL3
Protein Names	Recommended name: Four and a half LIM domains protein 3 Short name= FHL-3 Alternative name(s): Skeletal muscle LIM-protein 2 Short name= SLIM-2
Expression Region	1-280aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	58.1kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.