



# Recombinant Human Transcription factor 7 (TCF7)

<b>Product Code</b>	CSB-EP023305HU-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P36402
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPQLDSGGGG AGGGDDLGAPELLAFQDEG EEQDDKSRDS AAGPERDLAE LKSSLVNESE GAAGGAGIPG VPGAGAGARG EAEALGREHA AQRLFDPKLP EPLEDGLKAP ECTSGMYKET VYSAFNLLMH YPPPSGAGQH PQPQPPLHKA NQPPHGVPQL SLYEHFNPH PTPAPADISQ KQVHRPLQTP DLSGFYSLTS GSMGQLPHTV SWFTHPSLML GSGVPGHPAA IPHPAIVPPS GKQELQPFDR NLKTQAESKA EKEAKKPTIK KPLNAFMLYM KEMRAKVIAE CTLKESAAIN QILGRRWHAL SREEQAKYYE LARKERQLHM QLYPGWSARD NYGKKKRRSR EKHQESTTET NWPRELK DGN GQESLSMSSS SSPA
<b>Source</b>	E.coli
<b>Target Names</b>	TCF7
<b>Protein Names</b>	Recommended name: Transcription factor 7 Short name= TCF-7 Alternative name(s): T-cell-specific transcription factor 1 Short name= T-cell factor 1 Short name= TCF-1
<b>Expression Region</b>	1-384
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.