



# Recombinant Human Growth arrest-specific protein 6 (GAS6), Escherichia Coli

<b>Product Code</b>	CSB-EP619872HU
<b>Abbreviation</b>	GAS6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q14393
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>ALLPAREATQFLRPRQRRAFQVFEEAKQGHLERECVEELCSREEAREVFEND          PETDYFYPRYLDCINKYGSPYTKNSGFATCVQNLDPQCTPNPCDRKGTQACQ          DLMGNFFCLCKAGWGGRLCDKDVNECSQENGGCLQICHNKPGSFHCSCHS          GFELSSDGRTCQDIDECADSEACGEARCKNLPGSYSCLCDEGFAYSSQEKAC          RDVDECLQGRCEQVCVNSPGSYTCHCDGRGGLKLSQDMDTCEDILPCVPFS          VAKSVKSLYLGRMFSGTPVIRLRFKRLQPTLVAEFDRTFDPEGILLFAGGHQ          DSTWIVLALRAGRLELQLRYNGVGRVTSSGPVINHGMWQTISVEELARNLVIKV          NRDAVMKIAVAGDLFQPERGLYHLNLTVGGIPFHEKDLVQPINPRLDGCMRSW          NWLNGEDTTIQETVKVNTRMQCFVTERGSFYPGSGFAFYSLDYMRTPLDVG          TESTWEVEVVAHIRPAADTGVLFWAPDLRAVPLSVALVDYHSTKKLKKQLV          VLAVEHTALALMEIKVCDGQEHVVTVSLRDGEATLEVDGTRGQSEVSAAQLQE          RLAVLERHLRSPVLTFFAGGLPDVPVTSAPVTAFYRGCMTLEVNRRLLDLDEAA          YKHS DITAHSCPPVEPAAA</p>
<b>Source</b>	E.coli
<b>Target Names</b>	GAS6
<b>Protein Names</b>	Recommended name: Growth arrest-specific protein 6 Short name= GAS-6 Alternative name(s): AXL receptor tyrosine kinase ligand
<b>Expression Region</b>	31-678
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	<p>N-terminal His-tagged and C-terminal Myc-tagged          N-terminal His-SUMO-tagged and C-terminal Myc-tagged          N-terminal His-B2M-tagged and C-terminal Myc-tagged          N-terminal His-GST-tagged and C-terminal Myc-tagged          N-terminal His-tagged          N-terminal His-SUMO-tagged          N-terminal His-B2M-tagged          N-terminal His-GST-tagged</p>



N-terminal GST-tagged  
Tag-Free

---

**Protein Length**

Full Length of Mature Protein

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.