



Recombinant Human Growth arrest-specific protein 6 (GAS6), Escherichia Coli

Product Code	CSB-EP619872HU
Abbreviation	GAS6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14393
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	ALLPAREATQFLRPRQRRAFQVFEEAKQGHLERECVEELCSREEAREVFEND PETDYFYPRYLDKINKYKSPYTKNSGFATCVQNLDPDQCTPNPCDRKGTQACQ DLMGNFFCLCKAGWGGRLCDKDVNECSQENGGCLQICHNKPGSFHCSCHS GFELSSDGRTCQDIDECADSEACGEARCKNLPGSYSCLCDEGFAYSSQEKAC RDVDECLQGRCEQVCVNSPGSYTCHCDGRGGLKLSQDMDTCEDILPCVPFS VAKSVKSLYLGRMFSGTPVIRLRFKRLQPTLVAEFDRTFDPEGILLFAGGHQ DSTWIVLALRAGRLELQLRYNGVGRVTSSGPVINHGMWQTISVEELARNLVIKV NRDAVMKIAVAGDLFQPERGLYHLNLTVGIPFHEKDLVQPINPRLDGCMRSW NWLNGEDTTIQETVKVNTRMQCFVTERGSFYPGSGFAFYSLDYMRTPLDVG TESTWEVEVVAHIRPAADTGVLFALWAPDLRAVPLSVALVDYHSTKCLKKQLV VLAVEHTALALMEIKVCDGQEHVTVSLRDGEATLEVDGTRGQSEVSAAQLQE RLAVLERHLRSPVLTAFAGLPDVPVTSAPVTAFYRGCMTLEVNRRLLDLDEAA YKHS DITAHSCPPVEPAAA
Source	E.coli
Target Names	GAS6
Protein Names	Recommended name: Growth arrest-specific protein 6 Short name= GAS-6 Alternative name(s): AXL receptor tyrosine kinase ligand
Expression Region	31-678
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-SUMO-tagged and C-terminal Myc-tagged N-terminal His-B2M-tagged and C-terminal Myc-tagged N-terminal His-GST-tagged and C-terminal Myc-tagged N-terminal His-tagged N-terminal His-SUMO-tagged N-terminal His-B2M-tagged N-terminal His-GST-tagged



N-terminal GST-tagged
Tag-Free

Protein Length

Full Length of Mature Protein

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.