



# Recombinant Human Growth arrest-specific protein 6 (GAS6),Baculovirus

<b>Product Code</b>	CSB-BP619872HU
<b>Abbreviation</b>	GAS6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q14393
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ALLPAREATQFLRPRQRRAFQVFEEAKQGHLERECVEELCSREEAREVFEND PETDYFYPRYLDCINKYGSPYTKNSGFATCVQNLDPQCTPNPCDRKGTQACQ DLMGNFFCLCKAGWGGRLCDKDVNECSQENGGCLQICHNKPGSFHCSCHS GFELSSDGRTCQDIDECADSEACGEARCKNLPGSYSCLCDEGFAYSSQEKAC RDVDECLQGRCEQVCVNSPGSYTCHCDGRGGLKLSQDMDCEDILPCVPFS VAKSVKSLYLGRMFSGTPVIRLRFKRLQPTLVAEFDRTFDPEGILLFAGGHQ DSTWIVLALRAGRLELQLRYNGVGRVTSSGPVINHGMWQTISVEELARNLVIK NRDAVMKIAVAGDLFQPERGLYHLNLTVGGIPFHEKDLVQPINPRLDGCMRSW NWLNGEDTTIQETVKVNTRMQCFSVTERGSFYPGSGFAFYSLDYMRTPLDVG TESTWEVEVVAHIRPAADTGVLFWAPDLRAVPLSVALVDYHSTKCLKKQLV VLAVEHTALALMEIKVCDGQEHVVTVSLRDGEATLEVDGTRGQSEVSAAQLQE RLAVLERHLRSPVLTAFAGGLPDVPVTSAPVTAFYRGCMTLEVNRRLLDLDEAA YKHS DITAHSCPPVEPAAA
<b>Source</b>	Baculovirus
<b>Target Names</b>	GAS6
<b>Protein Names</b>	Recommended name: Growth arrest-specific protein 6 Short name= GAS-6Alternative name(s): AXL receptor tyrosine kinase ligand
<b>Expression Region</b>	31-678
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal His-tagged and C-terminal Myc-tagged N-terminal His-tagged Tag-Free
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.