




Recombinant Human Ficolin-2 (FCN2)

Product Code	CSB-EP620881HU
Abbreviation	Recombinant Human FCN2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15485
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	LQAADTCPEVKMVGLEGSDKLTILRGCPGLPGAPGPKGEAGTNGKRGERGPP GPPGKAGPPGPNAGPEPQPCLTGPRTCKDLLDRGHFLSGWHTIYLPDCRPL TVLCDMDTDGGGWTVFQRRVDGSDVDFYRDWATYKQGFGSRLGEFWLGNDN IHALTAQGTSELRVDLVDFEDNYQFAKYRSFKVADEAEKYNLVLGAFVEGSAG DSLTFHNNQSFSTKDQDNDLNTGNCAVMFQGAWWYKNCHVSNLNGRYLRG THGSFANGINWKSGKGYNYSYKVSEMKVRPA
Research Area	Cardiovascular
Source	E.coli
Target Names	FCN2
Protein Names	Recommended name: Ficolin-2 Alternative name(s): 37 kDa elastin-binding protein Collagen/fibrinogen domain-containing protein 2 EBP-37 Ficolin-B Ficolin-beta Hucolin L-ficolin Serum lectin p35
Expression Region	26-313aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	35.4 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.