



Recombinant Human Zinc finger protein 76 (ZNF76)

Product Code	CSB-EP027017HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36508
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MESLGLHTVT LSDGTTAYVQ QAVKGEKLL E GQVIQLEDGT TAYIHQVTVQ KEALSFEDGQ PVQLEDGSMA YIHRTPREGY DPSTLEAVQL EDGSTAYIHH PVAVPSESTI LAVQTEVGL E DLAAEDDEGF SADAVVALEQ YASKVLHDSQ IPRNGKGQV GDRAFRCGYK GCGRLYTTAH HMKVHERAHT GDRPYRCDFP SCGKAFATGY GLKSHVRTHT GEKPYKCPEE LCSKAFKTS G DLQKHVRTHT GERPFQCPFE GCGRSFTTSN IRKVHVRTHT GERPYTCPEP HCGRGFTSAT NYKNHVRIHT GEKPYVCTVP GCGKRFTEYS SLYKHHVVHT HCKPYTCSTC GKTYRQTSTL AMHKRSAHGE LEATEESEQA LYEQQQLEAA SAAEESPPPK RPRIAYLSEV KEERDDIPAQ VAMVTEEDGA PQVALITQDG AQQVSLSPED LQALGSAISM VTQHGSTTLT IPSPDADLAT SGTHTVTMVS ADGTQTQPVT IITSGAVVAE DSSVASLRHQ QVALLATANG THIAVQLEEQ QTLEEAINVA TAAMQQGAVT LETTVSESGC
Source	E.coli
Target Names	ZNF76
Protein Names	Recommended name: Zinc finger protein 76 Alternative name(s): Zinc finger protein 523
Expression Region	1-570
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.