



# Recombinant Human Zinc finger protein 76 (ZNF76)

<b>Product Code</b>	CSB-MP027017HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P36508
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MESLGLHTVT LSDGTTAYVQ QAVKGEKLL E GQVIQLEDGT TAYI HQVTVQ KEALSFEDGQ PVQLEDGSMA YIHRTPREGY DPSTLEAVQL EDGSTAYIHH PVAVPSESTI LAVQTEVGL E DLAAEDDEGF SADAVVALEQ YASKVLHDSQ IPRNGKGQQV GDRAFRCGYK GCGRLYTTAH HLKVHERAHT GDRPYRCDFP SCGKAFATGY GLKSHVRTHT GEKPYKCPEE LCSKAFKTS G DLQKHVRTHT GERPFQCPFE GCGRSFTTSN IRKVHVRTHT GERPYTCPEP HCGRGFTSAT NYKNHVRIHT GEKPYVCTVP GCGKRFTEYS SLYKHHVVHT HCKPYTCSTC GKTYRQTSTL AMHKRSAHGE LEATEESEQA LYEQQQLEAA SAAEESPPPK RPRIAYLSEV KEERDDIPAQ VAMVTEEDGA PQVALITQDG AQQVSLSPED LQALGSAISM VTQHGSTTLT IPSPDADLAT SGTHTVTMVS ADGTQTQPVT IITSGAVVAE DSSVASLRHQ QVALLATANG THIAVQLEEQ QTLEEAINVA TAAMQQGAVT LETTVSESGC
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ZNF76
<b>Protein Names</b>	Recommended name: Zinc finger protein 76 Alternative name(s): Zinc finger protein 523
<b>Expression Region</b>	1-570
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.