



# Recombinant Bovine Glutathione peroxidase 3 (GPX3)

<b>Product Code</b>	CSB-MP009868BO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P37141
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	QEKSKT DCHAGVGGTI YEYGALTIDG EEYIPFKQYA GKYILFVNVA SYUGLTGQYV ELNALQEELE PFGLVILGFP CNQFGKQEPG ENSEILATLK YVRPGGGFTP NFQLFEKGDV NGEKEQKFYT FLKNSCPPTS ELLGSPDRLF WPEMKVHDIR WNFEEKFLVGP DGIPIMRWYH RTTVNSVKMD ILTYMRRRAV WEAKGK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	GPX3
<b>Protein Names</b>	Recommended name: Glutathione peroxidase 3 Short name= GPx-3 Short name= GSHPx-3 EC= 1.11.1.9 Alternative name(s): Plasma glutathione peroxidase Short name= GPx-P Short name= GSHPx-P
<b>Expression Region</b>	25-226
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	This gene product belongs to the glutathione peroxidase family, which functions in the detoxification of hydrogen peroxide. It contains a selenocysteine (Sec) residue at its active site. The selenocysteine is encoded by the UGA codon, which normally signals translation termination. The 3' UTR of Sec-containing genes have a common stem-loop structure, the sec insertion sequence (SECIS), which is necessary for the recognition of UGA as a Sec codon rather than as a stop signal.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.