



Recombinant Escherichia coli O8 Aminomethyltransferase (gcvT)

Product Code	CSB-MP001684E00
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B7LYG9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O8 (strain IA11)
Purity	>85% (SDS-PAGE)
Sequence	MAQQTPLYEQ HTLCGARMVD FHGWMMPLHY GSQIDEHHAV RTDAGMFDVS HMTIVDLRGS RTREFRLRYLL ANDVAKLTKS GKALYSGMLN ASGGVIDDLI VYYFTEDFFR LVVNSATREK DLSWITQHAE PFGIEITVRD DLSMIAVQGP NAQAKAATLF NDAQRQAVEG MKPFFGVQAG DLFIATTGYT GEAGYEIALP NEKAADFWRALVEAGVKPCG LGARDTLRLE AGMNLYGQEM DETISPLAAN MGWTIAWEPA DRDFIGREAL EVQREHGTEK LVGLVMTEKG VLRNELPVRF TDAQGNQHEG IITSGTFSPT LGYSIALARV PEGIGETAIV QIRNREMPVK VTKPVFVRNG KAVA
Source	Mammalian cell
Target Names	gcvT
Protein Names	Recommended name: Aminomethyltransferase EC= 2.1.2.10 Alternative name(s): Glycine cleavage system T protein
Expression Region	1-364
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.