



Recombinant Escherichia coli O81 Hydroxyacylglutathione hydrolase (gloB)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP010121EOP |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | B7MQ21 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O81 (strain ED1a) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MNLNSIPAFD DNYIWVLNDE AGRCLIVDPG DAEPVLNAIS ANNWQPEAIF LTHHHHDHVG GVKELVEKFP QIVVYGPQET QDKGTTQVVK DGETAFVLGH EFSVIATPGH TLGHICYFSK PYLFCGDTLF SGGCGRLFEG TPSQMYQSLK KLSALPDDTL VCCAHEYTLS NMKFALSILP HDLSINDYYR KVKELRAKNQ ITLPVILKNE RQINVFLRTE NIDLINVINE ETLQPEER FAWLRSKKDR F |
| Source | E.coli |
| Target Names | gloB |
| Protein Names | Recommended name: Hydroxyacylglutathione hydrolase EC= 3.1.2.6 Alternative name(s): Glyoxalase II Short name= Glx II |
| Expression Region | 1-251 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |